

Male C57BL/6 Mouse Hepatocytes - Pooled, Cryopreserved

Product: M005152 | Lot: TNK

Storage: ≤ -150 °C

Post-Thaw Viability and Yield

Viability ≥ 75% post-thaw viability by trypan blue exclusion 76

8.63 million viable cells **Yield** ≥ 5 million viable cells

50 Donors Number of donors in pool

Post-thaw viability and yield are an average of samples from the beginning, middle, and end of the production run. All data generated without the use of Percoll.

Metabolic Activity

			Conc.		
	Enzyme	Substrate	[µM]	Metabolite	Result**
ECOD	7-Ethoxycoumarin O-deethylation	7-ethoxycoumarin	75	7-HC, 7-HCG, and 7-HCS*	199
UGT	7-Hydroxycoumarin glucuronidation	7-hydroxycoumarin	30	7-HCG	205
ST	7-Hydroxycoumarin sulfation	7-hydroxycoumarin	30	7-HCS	39.1

^{*7-}Hydroxycoumarin (7-HC), 7-Hydroxycoumarin glucuronide (7-HCG), 7-Hydroxycoumarin sulfate (7-HCS) **Metabolite rate of formation is measured in pmol/min/106 cells

Metabolic assays are run in triplicate. Activity results analyzed by HPLC-UV or LC/MS/MS validated procedures. Metabolite formation for all enzymes is measured after a 60 min. cell incubation at 37°C, 5% CO2 resulting in a final conc. of 1 million cells/mL.

Results for this lot have been derived through validated testing methods and confirmed by Quality Assurance.

Caution: This product is being sold for research and/or manufacturing purposes only. The biological samples supplied by BioIVT, or any material isolated from the samples, are for in-vitro research use only and are not to be used as a source of material for clinical therapies. Human material may be used in vivo in animals. The user assumes all responsibility for its usage and disposal, in accordance with all regulations.

M005152.03 | Page 1 of 1