

**Human Male Mongolia Pooled liver cytosol**

Product Number	LC-R-02M	Lot Number	QROX			
Gender	Male	Number of Animal in Pool	11			
Date of Preparation	2010-09	Storage Condition	-70°C			
<b>(1) Specification</b>						
Protein Content (mg/mL)			20			
Volume per vial (mL)			0.5			
Amount of protein per vial (mg/vial)			10			
<b>(2) In vitro intrinsic clearance (AOX)</b>						
Probe Substrate	( $\mu$ M)	$T_{1/2}$ (min)	Clint( $\mu$ L/min/mg protein)			
Vanillin	1.0	4.51	76.8			
<b>(3) DONOR INFORMATION</b>						
<b>Donor Number</b>	<b>Gender</b>	<b>Age</b>	<b>Race</b>	<b>Cause of Death</b>	<b>Smoker</b>	<b>Serology testing</b>
HM-03-006	Male	26	Mongolian	Trauma	Non smoker	Negative
HM-03-007	Male	24	Mongolian	Trauma	Non smoker	Negative
HM-03-008	Male	28	Mongolian	Trauma	Smoker	Negative
HM-03-011	Male	24	Mongolian	Trauma	Smoker	Negative
HM-03-012	Male	29	Mongolian	Trauma	Smoker	Negative
HM-03-013	Male	38	Mongolian	Trauma	Smoker	Negative
HM-03-014	Male	25	Mongolian	Trauma	Smoker	Negative
HM-03-015	Male	28	Mongolian	Trauma	Smoker	Negative
HM-02-001	Male	31	Mongolian	Trauma	Smoker	Negative
HM-02-002	Male	21	Mongolian	Trauma	Non smoker	Negative
HM-06-002	Male	35	Mongolian	Trauma	Smoker	Negative
Preparation by	Junhua Fan	Signature	<i>Junhua Fan</i>			
Reviewed by	Yi Sun	Signature	<i>Yi Sun</i>			
Approved by	Yuchang Mao	Signature	<i>yuchang mao</i>			