

## Porine Liver S9 – WZZS

Cat. No.: LS-XZ-02M

Lot No. : WZZS

Unit Size: 0.5 ml

Storage Condition: -70°C

Protein Concentration: 20 mg/ml

Updated at: 2011-04

### Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmol/mg/min)		
	Substrate ( $\mu$ M)	Mean	SD
Phenacetin deethylation	200	534	17.7
Testosterone 6 $\beta$ -hydroxylation	200	546	193

质量控制

李韶娟