## RILDResearch Institute for Liver Diseases (Shanghai) Co., Ltd.

## Male New Zealand Rabbit Liver S9

Product Number	LS-XBT-02	2M	Lot Number			ZXXZ
Gender	Male		Number of Animal in Pool		16	
Date of Preparation	2017-05		Storage Condition		-70°C	
(1) Specification						
Protein Content (mg/mL)					75	
Volume per vial (mL)					600	
Amount of protein per vial (g/bottle)					45	
(2) Enzyme Characterization						
Phenotyping Reaction		Substr	ate	Enzyme activity (pmole/mg/min)		
		(µM)	)	Mean		SD
Testosterone 6β-hydr	40		132		16.7	
7-hydroxylation CoumarinGlucuronidation		30		388		47.6
Prepared by	Weidong Yang	Signati	ure	your widong		
Reviewed by	Yi Sun	Signati	ure	Yi Sun Ferfer YU		
Approved by	Feifei Yu	Signati	ure	Ferfer YU		