

Beagle Dog Liver S9 – FBDN

Cat. No.: LS-BGQ-01M

Lot No.: FBDN

Unit Size: 0.5 mL

Storage Condition: -70 °C

Protein Concentration: 20 mg/ml

Updated at: 2011-4

Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmol/mg/min)		
	Substrate (μ M)	Mean	SD
Phenacetin deethylation	50	128	10.1
Testosterone 6 β -hydroxylation	100	70.8	8.83
Dextromethorphan O-demethylation	25	7.16	0.214
Chlorzoxazone 6-hydroxylation	200	166	44.7
Diclofenac 4-hydroxylation	100	16.7	0.882

Donor Information

Species	Beagle
Donor Number	1
Gender	Male

质量控制

李韶娟