

Wistar Rat Liver Cytosol- QYYH

Cat. No.: LC-DS-02M

Lot No.: QYYH

Unit Size: 0.5mL

Storage Condition: -70°C

Protein Concentration: 20 mg/ml

Updated: 2012-01

Lot Characterization Result

Phenotyping Reaction	Enzyme Activity (pmol/mg/min)		
	Substrate (μM)	Mean	SD
Phenacetin deethylation	50	0.513	0.378
Testosterone 6 β -hydroxylation	40	23.0	1.82

Donor Information

Species	Wistar rat
Donor Number	21
Gender	Male
Weight	220~250g

质量控制

李韶娟